



AT LISA toric 909M
AT LISA toric 909MP preloaded
Technical Specifications



We make it visible.

AT LISA toric 909M

AT LISA toric 909MP preloaded

AT LISA® toric 909M

Optic Design	Multifocal, diffractive, toric, +3.75 D near add at the IOL plane, aspheric (aberration correcting)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.5 mm
Company Labeled A-Constant ¹	118.3
Diopter Range	From -10.0 to +32.0 D Sphere: -10.0 to +32.0 D, 0.5 D increments Cylinder: +1.0 to +12.0 D, 0.5 D increments
ACD	5.14
Implantation in	Bag
Injector / Cartridge Set ²	For IOLs from -10.0 to +24.0 sphere with +1.0 to +4.0 cylinder*: AT.Shooter A2-2000 / ACM2 (1.5 mm) or VISCOJECT™ 1.8 Injector Set or Single-use injector A6 / *AT.Smart Cartridge Set (1.8 mm) * For diopter range outside -10.0 to +24.0 D sphere with +1.0 to +4.0 D cylinder, the correct injector / cartridge set will be provided with the IOL

AT LISA® toric 909MP preloaded

Optic Design	Multifocal, diffractive, toric, +3.75 D near add at the IOL plane, aspheric (aberration correcting)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	118.3
Diopter Range	From +6.0 to +24.0 D Sphere: +6.0 to +24.0 D, 0.5 D increments Cylinder: +1.0 to +4.0 D, 0.5 D increments
ACD	5.14
Implantation in	Bag
Injector / Cartridge Set ²	BLUEMIXS® 180

¹Please refer to our web pages for optimized A-Constants. ²Please refer to our web pages for the most up-to-date references.



0297 AT LISA 909M / 909MP
AT.Smart Cartridge Set
ACM2 Cartridge Set



0459 BLUEMIXS 180 injector



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.meditec.zeiss.com/iol
www.meditec.zeiss.com/contacts



Carl Zeiss Meditec SAS
Avenue Paul Langevin, BP5
17053 La Rochelle Cedex 9
France
www.meditec.zeiss.com/iol
www.meditec.zeiss.com/contacts

I/D01092_02/12_GB

AT LISA and BLUEMIXS are registered trademarks of Carl Zeiss Meditec. VISCOJECT is a trademark of Medical (CE 0482). The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper. © 2012 by Carl Zeiss Meditec AG. All copyrights reserved.



000000-1990-507

02/2012